

Northern Blot - Hiroki Kuroda

Sometimes useful to check amount and size of RNA

(RNA extraction with RNA-STAT)

- Mix with 800 μ l of RNA-Stat.
Note: 10 animal caps or five ventral marginal zone explants provide you almost same amount of total RNA as one whole embryo.
Note: You can keep sample in -80°C forever.
- Add 0.2 vol chloroform and vortex.
- Centrifuge at maximum speed (e.g. 15k rpm) for 10 min at 4°C .
- Recover the supernatant (usually 200 μ l) and do another chloroform extraction.
- Save 180 μ l of supernatant and add 180 μ l of isopropanol.
- Mix well and incubate in dry ice until frozen or in -80°C for 30 min.
Note: You can store sample in -20°C forever.
- Centrifuge at maximum speed for 15 min at 4°C .
- Wash the pellet with 70% EtOH and dissolve the pellet into 15 μ l of water.
Note: Use 4 μ l for the following room temperature reaction.

(Loading Sample)

Components	20 μ l final	Master Mix. (x20)
Total RNA Soln	5 μ l	----
0.4 mg/ml EtBr	1 μ l	20 μ l
20 x MOPS	0.5 μ l	10 μ l
formaldehyde	3.5 μ l	70 μ l
formamide	10 μ l	200 μ l

- Mix 15 μ l of master mix and 5 μ l of total RNA solution.
- Store at 60°C for 30 min.
- Store sample on ice.
- Add 2 μ l of Nuclease free 10 x Loading buffer.
- Electrophoresis using a formaldehyde gel.

(Blotting)

- Transfer overnight using the following system (From top to bottom)
 - Weight ($5\text{g}/\text{cm}^2$)
 - Glass plate
 - 5 cm of Kim Towel
 - 3 wet paper (11 x 6 cm)
 - wet Hybrid N+
 - Gel (gel surface-bottom)
 - 3 wet paper (11 x 20 cm)
 - glass, 20 x SSC
- Write a origin loading point on blotted paper.
- Check gel and membrane using UV to confirm whether blotting process is fine.
- Wash membrane with 2 x SSC for 1 min.
- UV crosslink and dry.

(Probe)

- Make DNA probe by PCR with α - ^{32}P -dCTP
- Purify with a PCR purification kit.

(20 x MOPS) 250 ml:

20.9 g of MOPS, 3.4 g $\text{CH}_3\text{COONa}(3\text{H}_2\text{O})$, 1.86 g of $\text{EDTA}(2\text{Na}/2\text{H}_2\text{O})$, pH 7.0

(20 x SSC) 500 ml

87.6 g of NaCl, 50 g of Sodium Citrate ($2\text{H}_2\text{O}$), pH 7.0

(Gel for Northern Blot)

-Mix 1 g of agarose, 5 ml of 20 x MOPS, and 77 ml of water using microwave
-Store at 60°C for 30 min.
-Add 18 ml of Formaldehyde

(Hybridization and detection)

Note: Use Hybri-bag for hybridization.

1. Wash the membrane with 5 x SSC (Note. RNA-side-up).
2. Move membrane into hybri-bag.
3. Add 15 ml pre-hybridization solution, then place the hybri-bag in the hybridization oven and incubate with rotation 6 hr at 42°C for DNA probe.
4. Pipet the desired volume of probe into hybri-bag.
Note. Double-stranded probe was denatured by heating in a water bath for 10 min at 100°C, then transferred to ice.
5. Continue to incubate with rotation overnight at 42°C for DNA probe.
6. Wash membrane twice for 10 min with wash-buffer at room temperature.
7. Wash membrane twice for 15 min with wash-buffer at 65°C.
8. Remove final wash solution and rinse membrane in 5 x SSC at room temperature.
9. Cover membrane in UV-transparent plastic wrap.
10. Do not allow membrane to dry out if it is to be probed again.
11. Blot was exposed at -80°C using Kodak XAR film and x-ray intensifying screens.