

RT-PCR – Hiroki Kuroda/Hojoon Lee

Most Common Technique for detecting molecular markers

RNA extraction with RNA-STAT

1. Lyse six embryos by pipetting with 800 μ l of RNA-Stat60.
- Note: Amount of total RNA is 3-4 μ g/embryo and 300-400 ng/AC.
2. Add 200 μ l of chloroform and vortex well.
3. Centrifuge at maximum speed (e.g. 15k rpm) for 10 min at 4 $^{\circ}$ C.
4. Recover the supernatant and repeat chloroform extraction.
5. Recover 380 μ l of the supernatant and add 380 μ l of isopropanol.
6. Mix well and incubate in dry ice until frozen or in -80 $^{\circ}$ C for 30 min.
- Note: You can store the sample in -80 $^{\circ}$ C forever.
7. Centrifuge at maximum speed for 15 min at 4 $^{\circ}$ C.
8. Wash the pellet with 70% EtOH and dissolve pellet in 20 μ l of water.
- Note. Use 4 μ l for the Reverse transcription reaction.

(Loading Sample)

Components	20 μ l final	Master Mix.
RNA	4 μ l	---
Water	6 μ l	300 μ l
5 x RT buffer	6 μ l	300 μ l
d[N] ₆ (100 mM)	3 μ l	150 μ l
dNTP (5 mM)	3 μ l	150 μ l
DTT (0.1M)	3 μ l	150 μ l
RNA guard	0.6 μ l	--- μ l
BSA (1 mg/ml)	3 μ l	150 μ l
MMLV RT enzyme	1.5 μ l	---

If you have 8 samples----


1. Mix 24 x 8.5 μ l of Master Mix, 0.6 x 8.5 μ l of RNA guard, and 1.5 x 8.5 μ l of MMLV RT enzyme very well (Called Master Enzyme Mix).
2. Mix 4 μ l of RNA solution and 26 μ l of Master Enzyme Mix very well.
3. Incubate at 42 $^{\circ}$ C for 1 h.
4. Add 30 μ l of water.
5. Store at -20 $^{\circ}$ C (Use as soon as possible. Do not keep for more than 3 months).

(PCR)

Components	20.4 μ l final	Master Mix.
cDNA solution	3 μ l	---
10 x PCR buffer (with Mg)	2 μ l	480 μ l
dNTP (5 mM)	0.8 μ l	192 μ l
α - ³² P-dCTP	0.1 μ l	---
Fw-Rv-Mixed primers (each 20 mM)	2 μ l	---
Water	12.2 μ l	2928 μ l
Taq polymerase	0.3 μ l	---

Major markers:

Stage 9	
Wnt target	<i>Chordin</i> <i>Noggin</i> <i>Xnr3</i> <i>Siamois</i> <i>Twin</i> <i>Pintallavis</i> <i>ADMP</i>
Nodal target	<i>Cerberus</i> <i>Xnr1</i> <i>Xnr2</i> <i>Xnr5</i> <i>Xnr6</i> <i>Sox17a</i> <i>Mix1</i> <i>Milk</i> <i>Mixer</i>
BMP target	<i>Vent1</i> <i>Vent2</i> <i>Msx1</i> <i>Msx2</i> <i>ID</i>
BCNE center	<i>Chordin</i> <i>Noggin</i> <i>Xnr3</i> <i>Siamois</i> <i>Twin</i> <i>Pintallavis</i> <i>ADMP</i>
Nieuwkoop center	<i>Cerberus</i> <i>Xnr1</i> <i>Xnr2</i> <i>Xnr5</i> <i>Xnr6</i>
Stage 13-14	Good for AC & VMZ assay
Cement Gland	<i>XAG</i>
Anterior Neural	<i>Otx2</i> <i>Six3</i>
Pan-neural	<i>NCAM</i> <i>Sox2</i>
Dorsal Mesodermal	<i>MyoD</i> <i>Myf5</i> <i>Xbra</i>
Epidermal	<i>Cytokeratin</i>
Ventral Mesodermal	<i>Bmp4</i> <i>Sizzled</i> <i>Crossveinless-2</i> <i>Vent1</i> <i>Vent2</i> <i>Msx1</i>



In the case that you have 8 different embryo samples and want to check 6 marker genes for each sample.
(Please note that font color scheme from the figure above is used in the following protocol) Set up the following cold reaction)

0-1. Set up PCR machine and prepare six sets of 8-connected PCR tubes.

0-2. Put 3 µl of cDNA soln into PCR tube and keep them in 4 °C.

0-3. A master mix pool will be prepared for each marker (Six in this example.) The master mix is prepared in 1.5 ml tubes (Hence six 1.5 ml tubes). Put 2 x 8.8 µl of respective primer pair mix into 1.5 ml tubes (total six tubes).

RADIOACTIVITY STARTS HERE!!! BE CAREFUL!!!

1. Mix 15 x 9 x 6 µl of PCR Master Mix, 0.3 x 9 x 6 µl of Taq polymerase, and 0.1 x 9 x 6 µl of α -³²P-dCTP (called PCR Master Enzyme Soln).
2. Take 15.4 x 8.8 µl of PCR Master Enzyme Soln, add to 1.5 ml tubes with primers (from 0-3). Mix very well (called PCR Master Enzyme & Primer Soln).
3. Add 17.4 µl of PCR Master Enzyme & Primer Soln (from #2 above) into each PCR tube (e.g. 17.4 µl each of Marker 1 into the 8 tubes of different samples)
4. Add one drop of mineral oil into each PCR tube.
5. Close tubes and start PCR.
6. Run on a 5 % acrylamide gel after PCR.

Components	Volume
Water	33 ml
40% acrylamide	5 ml
10x TBE	3.3 ml
APS	500 µl
TEMED	30 µl

5% acrylamide gel

	<i>Msx2</i> <i>Wnt8</i>
Stage 18-22	Good for VMZ assay
Ventral	<i>Sizzled</i> <i>β-Globin</i>
Paraxial	<i>PAPC</i>
Neural crest	<i>Slug</i>
Axial	<i>Shh</i> <i>Pintallavis</i>
Floor plate	<i>F-Spondin</i>
Dorsal	<i>α-Actin</i>
Cement Gland	<i>XAG</i>
Heart	<i>Nkx2.5</i>
Stage 22-26	Good for AC assay
Pan-neural	<i>NCAM</i> <i>N-Tubulin</i>
Anterior Neural	<i>Six3</i> <i>Otx2</i> <i>Pax6</i> <i>RX2a</i>
Mid-Hind-brain boundary	<i>Engrailed</i>
Hindbrain	<i>Krox20</i>
Posterior Neural	<i>Xlhbox6</i> (= <i>HoxB9</i>)
Mesoderm	<i>α-Actin</i> <i>β-Globin</i>
Cement Gland	<i>XAG</i>

Loading Control:

ODC: Always a good loading control marker

EF1α: Much stronger than *ODC* but not good before gastrula.