


## Spemann's organizer transplantation - Hiroki Kuroda

**Note: All surgery steps are done with forceps.**

**Note: Use 50 mm sized plastic plate.**

**Note: Use an embryo between stage 10-10.5 if you really want to have a complete secondary axis.**

- Put donor and host embryos together into 1x Steinberg's solution on a plate.
- Remove the chorion membrane from the dorsal animal side for host and ventral animal side for donor embryos.  
Note: If the embryos are significantly damaged, discard.
- Cut dorsal cap and remove yolk rich cells and the ectodermal region.
- Make slit at the ventral side of the host embryo. Note. The slit should reach to the blastocoel (this greatly helps healing).
- Transplant donor explant into ventral region in the host embryo.
- Leave transplanted embryos for 30 min (or until healed).
- Move transplanted embryo into a hole in 3% agarose plate in 0.1 x Steinberg's solution. Note: Do not use a high salt buffer, or you can get exogastrulation. Note: The hole in the agarose can be made with a heated glass ball at the tip of a glass pipet.
- Culture the embryo.



### 10 x Steinberg's solution

NaCl 34g (580 mM), KCl 0.5g (6.7 mM),  $\text{CaNO}_3 \cdot 4\text{H}_2\text{O}$  0.8g (3.4 mM),  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$  2g (8.3 mM), Kanamycin 0.1 g, Tris (MW=121) 6.0g (50 mM), adjust pH between 7.35-7.45 by HCl, adjust volume to 1L with  $\text{H}_2\text{O}$  and autoclave.

### 1 x Steinberg's solution

Kanamycin 0.1g, 10x Steinberg's solution 100 ml, adjust volume to 1L  $\text{H}_2\text{O}$  and autoclave.